

## anti- TRIM31 antibody

### Product Information

Catalog No.:	FNab08979
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Regulator of Src-induced anchorage independent cell growth(By similarity). May have E3 ubiquitin-protein ligase activity.

### Immunogen information

Immunogen:	tripartite motif-containing 31
Synonyms:	C6orf13, HCG1, HCGI, RNF, TRIM31, tripartite motif containing 31
Observed MW:	48 kDa
UniprotID :	Q9BZY9

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human colon cancer using FNab08979(TRIM31 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab08979(TRIM31 antibody) at dilution of 1:400