

## anti- TRPM1 antibody

### Product Information

Catalog No.:	FNab09020
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Cation channel essential for the depolarizing photoresponse of retinal ON bipolar cells. It is part of the GRM6 signaling cascade. May play a role in metastasis suppression(By similarity). May act as a spontaneously active, calcium-permeable plasma membrane channel.

### Immunogen information

Immunogen:	transient receptor potential cation channel, subfamily M, member 1
Synonyms:	LTRPC1, Melastatin 1, MLSN, MLSN1, TRPM1
Observed MW:	182 kDa
Uniprot ID :	Q7Z4N2

### Application

Reactivity:	Human
Tested Application:	ELISA, IF
Recommended dilution:IF:	1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)