

anti- TRPM6 antibody

Product Information

Catalog No.:	FNab09022
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Essential ion channel and serine/threonine-protein kinase. Crucial for magnesium homeostasis. Has an important role in epithelial magnesium transport and in the active magnesium absorption in the gut and kidney. Isoforms of the type M6-kinase lack the ion channel region.

Immunogen information

Immunogen:	transient receptor potential cation channel, subfamily M, member 6
Synonyms:	CHAK2, Channel kinase 2, HMGX, HOMG, HOMG1, HSH, TRPM6
Observed MW:	232 kDa
Uniprot ID :	Q9BX84

Application

Reactivity:	Human
Tested Application:	ELISA, IF
Recommended dilution:IF:	1:10-1:100
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com