

anti- TFSM antibody

Product Information

Catalog No.:	FNab09039
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Associates with the EF-Tu.GDP complex and induces the exchange of GDP to GTP. It remains bound to the aminoacyl-tRNA.EF-Tu.GTP complex up to the GTP hydrolysis stage on the ribosome.

Immunogen information

Immunogen:	Ts translation elongation factor, mitochondrial
Synonyms:	COXPD3, EF TS, EF Tsmt, TFSM
Observed MW:	31 kDa, Observed 39 kDa
UniprotID :	P43897

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon cancer using FNab09039(TSFM antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09039(TSFM Antibody) at dilution of 1:600