

## anti- TSPAN15 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab09058   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

None

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | tetraspanin 15  |
| Synonyms:    | 2700063A19Rik, NET 7, NET7, Tetraspan NET 7, tetraspanin 15, TM4SF15, Tspan 15, TSPAN15 |
| Observed MW: | 33 kDa  |
| UniprotID :  | O95858  |

### Application


|                       |                                   |
|-----------------------|-----------------------------------|
| Reactivity:           | Human                             |
| Tested Application:   | ELISA, WB, IHC                    |
| Recommended dilution: | WB: 1:200-1:2000; IHC: 1:20-1:200 |
| Image:                |                                   |

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab09058(TSPAN15 antibody) at dilution of 1:100



human ovary cancer were subjected to SDS PAGE followed by western blot with FNab09058(TSPAN15 antibody) at dilution of 1:500