

## anti- TSPAN3 antibody

### Product Information

Catalog No.:	FNab09060
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Regulates the proliferation and migration of oligodendrocytes, a process essential for normal myelination and repair.

### Immunogen information

Immunogen:	tetraspanin 3
Synonyms:	Tetraspanin 3, Tetraspanin TM4 A, TM4 A, TM4SF8, TSPAN 3, TSPAN3
Observed MW:	28 kDa
Uniprot ID :	O60637

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)