

## anti- TSPY3 antibody

### Product Information

Catalog No.:	FNab09067
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Unknown. May be involved in sperm differentiation and proliferation(By similarity).

### Immunogen information

Immunogen:	testis specific protein, Y-linked 3
Synonyms:	Cancer/testis antigen 78, CT78, TSPY, TSPY1, TSPY3
Observed MW:	35 kDa, Observed 40 kDa
UniprotID :	P0CV98

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:1000-1:10000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis using FNab09067(TSPY3 antibody) at dilution of 1:100



LNCaP cells were subjected to SDS PAGE followed by western blot with FNab09067(TSPY3 antibody) at dilution of 1:2000