

anti- TTC30B antibody

Product Information

Catalog No.:	FNab09080
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Required for polyglutamylation of axonemal tubulin. Plays a role in anterograde intraflagellartransport(IFT), the process by which cilia precursors are transported from the base of the cilium to the site of their incorporation at the tip.

Immunogen information

Immunogen:	tetratricopeptide repeat domain 30B
Synonyms:	FLJ30990, TPR repeat protein 30B, TTC30B
Observed MW:	57 kDa
UniprotID :	Q8N4P2

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09080(TTC30B Antibody) at dilution of 1:600