

anti- TXNL1 antibody

Product Information

Catalog No.:	FNab09130
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Immunogen information

Immunogen:	thioredoxin-like 1
Synonyms:	Thioredoxin like 1, Thioredoxin like protein 1, TRP32, Tx1, TXL 1, TXNL, TXNL1
Observed MW:	30 kDa
UniprotID :	O43396

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain using FNab09130(TXNL1 antibody) at dilution of 1:100



K-562 cells were subjected to SDS PAGE followed by western blot with FNab09130(TXNL1 antibody) at dilution of 1:1000