

anti- UAP1L1 antibody

Product Information

Catalog No.:	FNab09141
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	UDP-N-acteylglucosaminepyrophosphorylase 1-like 1
Synonyms:	UDP-N-acetylhexosaminepyrophosphorylase-like protein 1,UAP1L
Observed MW:	57 kDa
UniprotID :	Q3KQV9

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:WB:	1:200-1:2000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with FNab09141(UAP1L1 Antibody) at dilution of 1:600