

anti- UBE2U antibody

Product Information

Catalog No.:	FNab09186
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Immunogen information

Immunogen:	ubiquitin-conjugating enzyme E2U (putative)
Synonyms:	RP4 636O23.1, UBE2U, Ubiquitin carrier protein U, Ubiquitin protein ligase U
Observed MW:	42 kDa
UniprotID :	Q5VVX9

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver using FNab09186(UBE2U antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab09186(UBE2U antibody) at dilution of 1:1000