

anti- UBQLNL antibody

Product Information

Catalog No.:	FNab09208
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ubiquilin-like
Synonyms:	Ubiquilin like, Ubiquilin like protein, UBQLNL
Observed MW:	53 kDa, Observed 70 kDa
UniprotID :	Q8IYU4

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1: 10-1:100; IHC: 1:20-1:200; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab09208(UBQLNL Antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09208(UBQLNL antibody) at dilution of 1:500