

anti- UHRF1BP1L antibody

Product Information

Catalog No.:	FNab09246
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	UHRF1 binding protein 1-like
Synonyms:	KIAA0701, UHRF1 binding protein 1 like, UHRF1BP1L
Observed MW:	164 kDa
Uniprot ID :	A0JNW5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:500-1:5000; IP: 1: 500-1: 5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP Result of anti-UHRF1BP1L (IP:FNab09246, 5ug; Detection:FNab09246 1:1000) with mouse brain tissue lysate 5000ug.



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09246(UHRF1BP1L Antibody) at dilution of 1:1500