

anti- UNC45B antibody

Product Information

Catalog No.:	FNab09259
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as a co-chaperone for HSP90 and is required for proper folding of the myosin motor domain. Plays a role in sarcomere formation during muscle cell development. Is necessary for normal early lens development.

Immunogen information

Immunogen:	unc-45 homolog B(C. elegans)
Synonyms:	CMYA4, FLJ38610, Protein unc 45 homolog B, SMUNC45, unc 45 homolog B(C. elegans), Unc 45B, UNC45, UNC45B
Observed MW:	95 kDa
Uniprot ID :	Q8IWX7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab09259(UNC45B antibody) at dilution of 1:500