

## anti- UPK3B antibody

### Product Information

|               |                                                                                                                 |
|---------------|-----------------------------------------------------------------------------------------------------------------|
| Catalog No.:  | FNab09276                                                                                                       |
| Size:         | 100µg                                                                                                           |
| Form:         | liquid                                                                                                          |
| Purification: | Immunogen affinity purified                                                                                     |
| Purity:       | ≥95% as determined by SDS-PAGE                                                                                  |
| Host:         | Rabbit                                                                                                          |
| Clonality:    | polyclonal                                                                                                      |
| Clone ID:     | None                                                                                                            |
| IsoType:      | IgG                                                                                                             |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

Component of the asymmetric unit membrane(AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in AUM-cytoskeleton interaction in terminally differentiated urothelial cells. It also contributes to the formation of urothelial glycocalyx which may play an important role in preventing bacterial adherence(By similarity).

### Immunogen information

|              |                                                                |
|--------------|----------------------------------------------------------------|
| Immunogen:   | urolakin 3B                                                    |
| Synonyms:    | FLJ32198, P35, UP3b, UPIIIb, UPK3B, urolakin 3B, UroplakinIIIb |
| Observed MW: | 34 kDa                                                         |
| UniprotID :  | Q9BT76                                                         |

### Application

|                       |                                   |
|-----------------------|-----------------------------------|
| Reactivity:           | Human, Mouse                      |
| Tested Application:   | ELISA, IHC, IP                    |
| Recommended dilution: | IHC: 1:20-1:200; IP: 1:200-1:1000 |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Image:



IP Result of anti-UPK3B (IP:FNab09276, 2ug;  
Detection:FNab09276 1:300) with mouse testis  
tissue lysate 4000ug.

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)