

anti- USP24 antibody

Product Information

Catalog No.:	FNab09320
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Deubiquitinates DDB2, preventing its proteasomal degradation.

Immunogen information

Immunogen:	ubiquitin specific peptidase 24
Synonyms:	Deubiquitinating enzyme 24, FLJ31309, KIAA1057, Ubiquitin thiolesterase 24, USP24
Observed MW:	290 kDa
Uniprot ID :	Q9UPU5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab09320(USP24 antibody) at dilution of 1:300