

anti- USP32 antibody

Product Information

Catalog No.:	FNab09324
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ubiquitin specific peptidase 32
Synonyms:	Deubiquitinating enzyme 32, NY REN 60, Ubiquitin thiolesterase 32, USP10, USP32
Observed MW:	182 kDa
Uniprot ID :	Q8NFA0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IP: 1: 200-1: 2000; IHC: 1: 20-1: 200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas using FNab09324(USP32 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09324(USP32 antibody) at dilution of 1:400