

## anti- USP45 antibody

### Product Information

Catalog No.:	FNab09335
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	ubiquitin specific peptidase 45
Synonyms:	Deubiquitinating enzyme 45, Ubiquitin thiolesterase 45, USP45
Observed MW:	92 kDa, Observed 45 kDa
Uniprot ID :	Q70EL2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1: 50-1: 500
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab09335(USP45 antibody) at dilution of 1:1500