

## anti- UTP15 antibody

### Product Information

Catalog No.:	FNab09348
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Involved in nucleolar processing of pre-18S ribosomal RNA.

### Immunogen information

Immunogen:	UTP15, U3 small nucleolar ribonucleoprotein, homolog(S. cerevisiae)
Synonyms:	FLJ12787, NET21, UTP15
Observed MW:	58 kDa,37 kDa
Uniprot ID :	Q8TED0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab09348(UTP15 antibody) at dilution of 1:300