

anti- UTP3 antibody

Product Information

Catalog No.:	FNab09350
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain(By similarity).

Immunogen information

Immunogen:	UTP3, small subunit(SSU) processome component, homolog(S. cerevisiae)
Synonyms:	UTP3, FLJ23256, CRLZ1, SAS10, UTP3 homolog
Observed MW:	55 kDa, Observed 30 kDa, 70 kDa
UniprotID :	Q9NQZ2

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:10000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



A549 cells were subjected to SDS PAGE followed by western blot with FNab09350(UTP3 antibody) at dilution of 1:1000