

## anti- VARS2 antibody

### Product Information

Catalog No.:	FNab09372
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	valyl-tRNA synthetase 2, mitochondrial(putative)
Synonyms:	KIAA1885, Valine tRNA ligase, ValRS, Valyl tRNA synthetase like, VARS2, VARS2L, VARSL
Observed MW:	118 kDa
Uniprot ID :	Q5ST30

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:50-1:500; IP: 1:200-1:2000
Image:	

1


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab09372(VARS2 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09372(VARS2 antibody) at dilution of 1:500