

anti- VEZT antibody

Product Information

Catalog No.:	FNab09396
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a pivotal role in the establishment of adherens junctions and their maintenance in adult life. In case of Listeria infection, promotes bacterial internalization by participating in myosin VIIa recruitment to the entry site.

Immunogen information

Immunogen:	vezatin, adherens junctions transmembrane protein
Synonyms:	VEZT, VEZATIN, DKFZp761C241
Observed MW:	89 kDa, 125 kDa
Uniprot ID :	Q9HBM0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain using FNab09396(VEZT antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09396(VEZT antibody) at dilution of 1:500