

anti- VGLUT1 antibody

Product Information

Catalog No.:	FNab09400
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	solute carrier family 17(sodium-dependent inorganic phosphate cotransporter), member 7
Synonyms:	Brain-specific Na(+)-dependent inorganic phosphate cotransporter,Solute carrier family 17 member 7,BNPI, VGLUT1,SLC17A7
Observed MW:	62 kDa
Uniprot ID :	Q9P2U7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, FC
Recommended dilution:	FC:N/A; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemical of paraffin-embedded human brain using FNab09400(SLC17A7 antibody) at dilution of 1:50