
anti- VPR antibody

Product Information

Catalog No.:	FNab09423
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

This antibody recognize the VPR gene of Human immunodeficiency virus 1.

Immunogen information

Immunogen:	VPR
Synonyms:	Protein Vpr, Viral protein R, R ORF protein
Uniprot ID :	P05928

Application

Reactivity:	Human immunodeficiency virus 1
Tested Application:	ELISA

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com