

## anti- VPREB3 antibody

### Product Information

Catalog No.:	FNab09425
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Associates with the Ig-mu chain to form a molecular complex that is expressed on the surface of pre-B-cells.

### Immunogen information

Immunogen:	pre-B lymphocyte 3
Synonyms:	8HS20, N27C7 2, pre B lymphocyte 3, Pre B lymphocyte protein 3, Protein VPreB3, VPREB3
Observed MW:	14 kDa, observed 46 kDa
UniprotID :	Q9UKI3

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab09425(VPREB3 antibody) at dilution of 1:400