

## anti- VPS26A antibody

### Product Information

Catalog No.:	FNab09432
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Acts as component of the retromer cargo-selective complex(CSC). The CSC is believed to be the core functional component of retromer or respective retromer complex variants acting to prevent missorting of selected transmembrane cargo proteins into the lysosomal degradation pathway. The recruitment of the CSC to the endosomal membrane involves RAB7A and SNX3. The SNX-BAR retromer mediates retrograde transport of cargo proteins from endosomes to the trans-Golgi network(TGN) and is involved in endosome-to-plasma membrane transport for cargo protein recycling. The SNX3-retromer mediates the retrograde endosome-to-TGN transport of WLS distinct from the SNX-BAR retromer pathway. The SNX27-retromer is believed to be involved in endosome-to-plasma membrane trafficking and recycling of a broad spectrum of cargo proteins(Probable). The CSC seems to act as recruitment hub for other proteins, such as the WASH complex and TBC1D5(Probable). Required for retrograde transport of lysosomal enzyme receptor IGF2R(PubMed:15078902, PubMed:15078903). Required to regulate transcytosis of the polymeric immunoglobulin receptor(pIgR-pIgA)(PubMed:15247922). Required for the endosomal localization of FAM21A(indicative for the WASH complex)(PubMed:22070227). Required for the endosomal localization of TBC1D5(PubMed:20923837). Mediates retromer cargo recognition of SORL1 and is involved in trafficking of SORL1 implicated in sorting and processing of APP(PubMed:22279231). Involved in retromer-independent lysosomal sorting of F2R(PubMed:16407403). Involved in recycling of ADRB2(PubMed:21602791). Enhances the affinity of SNX27 for PDZ-binding motifs in cargo proteins(By similarity).

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)


---

**Immunogen information**

Immunogen: vacuolar protein sorting 26 homolog A  
Synonyms: FLJ12930, HB58, Hbeta58, hVPS26, PEP8A, Vesicle protein sorting 26A, VPS26, VPS26A  
Observed MW: 38 kDa  
Uniprot ID : O75436

**Application**

Reactivity: Human, Mouse, Rat  
Tested Application: ELISA, WB, IHC,IF  
Recommended dilution:WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200  
Image:



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab09432(VPS26A Antibody) at dilution of 1:200


**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09432(VPS26A Antibody) at dilution of 1:300