

## anti- WDR76 antibody

### Product Information

Catalog No.:	FNab09502
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Specifically binds 5-hydroxymethylcytosine(5hmC), suggesting that it acts as a specific reader of 5hmC.

### Immunogen information

Immunogen:	WD repeat domain 76
Synonyms:	CDW14, FLJ12973, WD repeat domain 76, WDR76
Observed MW:	65 kDa
UniprotID :	Q9H967

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09502(WDR76 Antibody) at dilution of 1:300