

anti- WNT3 antibody

Product Information

Catalog No.:	FNab09518
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube (By similarity).

Immunogen information

Immunogen:	wingless-type MMTV integration site family, member 3
Synonyms:	INT4, Proto oncogene Int 4 homolog, Proto oncogene Wnt 3, WNT3
Observed MW:	38 kDa
UniprotID :	P56703

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone, Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab09518(WNT3 antibody) at dilution of 1:600