

anti- WRNIP1 antibody

Product Information

Catalog No.:	FNab09525
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Functions as a modulator for initiation or reinitiation events during DNA polymerase delta-mediated DNA synthesis. Has an intrinsic ATPase activity that functions as a sensor of DNA damage or of arrested replication forks and regulates the extent of DNA synthesis.

Immunogen information

Immunogen:	Werner helicase interacting protein 1
Synonyms:	WRNIP1, WHIP, bA420G6.2, RP11 420G6.2, ATPase WRNIP1
Observed MW:	50 kDa
Uniprot ID :	Q96S55

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



human colon tissue were subjected to SDS PAGE followed by western blot with FNab09525(WRNIP1 antibody) at dilution of 1:500