

anti- WTAP antibody

Product Information

Catalog No.:	FNab09529
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Regulatory subunit of the WMM N6-methyltransferase complex, a multiprotein complex that mediates N6-methyladenosine(m6A) methylation of some adenosine residues of some mRNAs and plays a role in the efficiency of mRNA splicing, processing and mRNA stability. Required for accumulation of METTL3 and METTL14 to nuclear speckle(PubMed:24316715, PubMed:24407421, PubMed:24981863). Acts as a mRNA splicing regulator(PubMed:12444081). Regulates G2/M cell-cycle transition by binding to the 3' UTR of CCNA2, which enhances its stability(PubMed:17088532). Impairs WT1 DNA-binding ability and inhibits expression of WT1 target genes(PubMed:17095724).

Immunogen information

Immunogen:	Wilms tumor 1 associated protein
Synonyms:	DKFZp686F20131, Female lethal(2)D homolog, hFL(2)D, KIAA0105, WT1 associated protein, WTAP
Observed MW:	55kDa
Uniprot ID :	Q15007

Application

1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com

Reactivity: Human, Mouse, Rat


Tested Application: ELISA, IP, WB

Recommended dilution: WB: 1:500-1:2000

Image:



IP Result of anti-WTAP (IP:FNab09529, 3ug;
Detection:FNab09529 1:500) with HEK-293 cells
lysate 2500ug.



Jurkat cells were subjected to SDS PAGE followed
by western blot with FNab09529(WTAP antibody)
at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com