

## anti- WTX antibody

### Product Information

Catalog No.:	FNab09531
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	family with sequence similarity 123B
Synonyms:	AMER1, FAM123B, Protein FAM123B, WTX
Observed MW:	134 kDa
UniprotID :	Q5JTC6

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab09531(FAM123B Antibody) at dilution of 1:200