

anti- XPNPEP1 antibody

Product Information

Catalog No.:	FNab09548
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Contributes to the degradation of bradykinin. Catalyzes the removal of a penultimate prolyl residue from the N-termini of peptides, such as Arg-Pro-Pro.

Immunogen information

Immunogen:	X-prolyl aminopeptidase(aminopeptidase P) 1, soluble
Synonyms:	XPNPEP1, XPNPEPL, Xaa Pro aminopeptidase 1, Cytosolic aminopeptidase P, XPNPEP
Observed MW:	
Uniprot ID :	Q9NQW7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA
Recommended dilution:	None

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com