

anti- YPEL3 antibody

Product Information

Catalog No.:	FNab09569
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in proliferation and apoptosis in myeloid precursor cells.

Immunogen information

Immunogen:	yippee-like 3(Drosophila)
Synonyms:	Protein yippee like 3, yippee like 3(Drosophila), YPEL3
Observed MW:	14 kDa
Uniprot ID :	P61236

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab09569(YPEL3 antibody) at dilution of 1:1000