

anti- YY1AP1 antibody

Product Information

Catalog No.:	FNab09579
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Enhances transcription activation by YY1. May play a role in cell cycle regulation.

Immunogen information

Immunogen:	YY1 associated protein 1
Synonyms:	FLJ10875, FLJ13914, HCCA1, HCCA2, YAP, YY1 associated protein 1, YY1AP, YY1AP1
Observed MW:	60 kDa, 75 kDa
Uniprot ID :	Q9H869

Application

Reactivity:	Human
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



MDA-MB-231 cells were subjected to SDS PAGE followed by western blot with FNab09579(YY1AP1 antibody) at dilution of 1:800