

anti- ZC3H11A antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab09600 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Involved in nuclear mRNA export; probably mediated by association with the TREX complex.

Immunogen information

| | |
|--------------|---|
| Immunogen: | zinc finger CCCH-type containing 11A |
| Synonyms: | DKFZp686D03108, DKFZp686F14109, DKFZp781G2455, KIAA0663, ZC3H11A, ZC3HDC11A |
| Observed MW: | 89 kDa |
| Uniprot ID : | O75152 |

Application

| | |
|-----------------------|----------------------------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution: | WB: 1:200-1:2000; IF: 1:20-1:200 |
| Image: | |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab09600(ZC3H11A Antibody) at dilution of 1:600