

## anti- ZDHHC15 antibody

### Product Information

Catalog No.:	FNab09616
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Palmitoyltransferase specific for GAP43 and DLG4/PSD95.

### Immunogen information

Immunogen:	zinc finger, DHHC-type containing 15
Synonyms:	DHHC 15, FLJ31812, MRX91, Palmitoyltransferase ZDHHC15, ZDHHC15
Observed MW:	42 kDa
Uniprot ID :	Q96MV8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with FNab09616(ZDHHC15 Antibody) at dilution of 1:1000