

## anti- ZDHHC5 antibody

### Product Information

Catalog No.:	FNab09617
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Palmitoyl acyltransferase for the G-protein coupled receptor SSTR5. Also palmitoylates FLOT2 (By similarity).

### Immunogen information

Immunogen:	zinc finger, DHHC-type containing 5
Synonyms:	DHHC 5, DKFZp586K0524, KIAA1748, ZDHHC5, Zinc finger protein 375, ZNF375
Observed MW:	78 kDa
UniprotID :	Q9C0B5

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunofluorescent analysis of HeLa cells, using FNab09617(ZDHHC5 antibody) at 1:25 dilution.



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab09617(ZDHHC5 antibody) at dilution of 1:500