

anti- ZDHHC7 antibody

Product Information

Catalog No.:	FNab09618
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Palmitoyltransferase with broad specificity. Palmitoylates SNAP25 and DLG4/PSD95(By similarity). Palmitoylates sex steroid hormone receptors, including ESR1, PGR and AR, thereby regulating their targeting to the plasma membrane and their function in rapid intracellular signaling upon binding of sex hormones(PubMed:22031296). May play a role in follicle stimulation hormone(FSH) activation of testicular Sertoli cells(By similarity).

Immunogen information

Immunogen:	zinc finger, DHHC-type containing 7
Synonyms:	DHHC 7, FLJ10792, FLJ20279, Palmitoyltransferase ZDHHC7, ZDHHC7, Zinc finger protein 370, ZNF370
Observed MW:	45 kDa
Uniprot ID :	Q9NXF8

Application

Reactivity:	Human
Tested Application:	ELISA, WB

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Recommended dilution: WB: 1:500-1:5000

Image:



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09618(ZDHHC7 Antibody) at dilution of 1:1000