

anti- ZIM3 antibody

Product Information

Catalog No.:	FNab09638
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger, imprinted 3
Synonyms:	ZIM3, Zinc finger imprinted 3, zinc finger, imprinted 3, ZNF657
Observed MW:	60 kDa
Uniprot ID :	Q96PE6

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab09638(ZIM3 Antibody) at dilution of 1:300