

anti- ZKSCAN1 antibody

Product Information

Catalog No.:	FNab09645
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger with KRAB and SCAN domains 1
Synonyms:	9130423L19Rik, KOX18, PHZ 37, Zinc finger protein 139, Zinc finger protein 36, Zinc finger protein KOX18, ZKSCAN1, ZNF139, ZNF36
Observed MW:	62 kDa
UniprotID :	P17029

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09645(ZKSCAN1 Antibody) at dilution of 1:600