

anti- ZNF174 antibody

Product Information

Catalog No.:	FNab09661
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcriptional repressor.

Immunogen information

Immunogen:	zinc finger protein 174
Synonyms:	ZNF174, AW 1, zinc finger protein 174, ZSCAN8
Observed MW:	46 kDa
Uniprot ID :	Q15697

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IF: 1:10-1:100; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis using FNab09661(ZNF174 antibody) at dilution of 1:100



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab09661(ZNF174 antibody) at dilution of 1:500