

anti- ZNF271 antibody

Product Information

Catalog No.:	FNab09678
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May act to control gene activity during the pachytene stage of meiotic prophase. May function as a transcription activator.

Immunogen information

Immunogen:	zinc finger protein 271
Synonyms:	HZF7, zinc finger protein 271, ZNF271, ZNFEB
Observed MW:	48 kDa
Uniprot ID :	P15620

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab09678(ZNF271 antibody) at dilution of 1:1000