

anti- ZNF326 antibody

Product Information

Catalog No.:	FNab09686
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Core component of the DBIRD complex, a multiprotein complex that acts at the interface between core mRNP particles and RNA polymerase II(RNAPII) and integrates transcript elongation with the regulation of alternative splicing: the DBIRD complex affects local transcript elongation rates and alternative splicing of a large set of exons embedded in(A + T)-rich DNA regions. May play a role in neuronal differentiation and is able to bind DNA and activate expression in vitro.

Immunogen information

Immunogen:	zinc finger protein 326
Synonyms:	DJ871E2.1, FLJ20403, ZAN75, Zfp326, zinc finger protein 326, ZNF326
Observed MW:	67 kDa
UniprotID :	Q5BKZ1

Application

Reactivity:	Human, Mouse
-------------	--------------

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Tested Application: ELISA, WB, IF

Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:500; IF: 1:10-1:100

Image:



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09686(ZNF326 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com