

## anti- Znt5 antibody

### Product Information

Catalog No.:	FNab09748
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	solute carrier family 30(zinc transporter), member 5
Synonyms:	FLJ12496, FLJ12756, hZTL1, SLC30A5, Zinc transporter 5, ZnT 5, ZnT like transporter 1, ZNT5, ZNTL1, ZTL1
Observed MW:	84 kDa
UniprotID :	Q8TAD4

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-Znt5 (IP:FNab09748, 4ug;  
Detection:FNab09748 1:500) with mouse liver  
tissue lysate 4000ug.