

anti- ZPBP antibody

Product Information

Catalog No.:	FNab09754
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be implicated in gamete interaction during fertilization.

Immunogen information

Immunogen:	zona pellucida binding protein
Synonyms:	SP38, zona pellucida binding protein, ZPBP, ZPBP1
Observed MW:	40 kDa, observed 34 kDa
UniprotID :	Q9BS86

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:200-1:1000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-ZPBP (IP:FNab09754, 3ug;
Detection:FNab09754 1:300) with mouse testis
tissue lysate 2800ug.