

anti- ZRANB2 antibody

Product Information

Catalog No.:	FNab09755
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Splice factor required for alternative splicing of TRA2B/SFRS10 transcripts. May interfere with constitutive 5'-splice site selection.

Immunogen information

Immunogen:	zinc finger, RAN-binding domain containing 2
Synonyms:	DKFZp686J1831, DKFZp686N09117, FLJ41119, Zinc finger protein 265, Zinc finger, splicing, ZIS, ZIS1, ZIS2, ZNF265, ZRANB2
Observed MW:	50 kDa
Uniprot ID :	O95218

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab09755(ZRANB2 Antibody) at dilution of 1:600