

## anti- ZSCAN20 antibody

### Product Information

Catalog No.:	FNab09757
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May be involved in transcriptional regulation.

### Immunogen information

Immunogen:	zinc finger and SCAN domain containing 20
Synonyms:	KOX29, ZFP 31, Zinc finger protein 31, Zinc finger protein 360, Zinc finger protein KOX29, ZNF31, ZNF360, ZSCAN20
Observed MW:	110 kDa
Uniprot ID :	P17040

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab09757(ZSCAN20 Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09757(ZSCAN20 antibody) at dilution of 1:1500