
anti- Ferritin antibody

Product Information

Catalog No.:	FNab09810
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4A9
IsoType:	IgG1
Storage:	20mMPBS, pH7.4 with 0.1% proclin300, 2-8°C for 2 weeks, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Ferritin
Observed MW:	170 kDa
Uniprot ID :	P02794

Application

Reactivity:	Human
Application:	ELISA, CLIA, GICA (capture antibody)

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com