
JAK1 antibody

Product Information


Catalog No.:	FNab09903
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Janus kinase 1
Synonyms:	JAK 1, JAK1, JAK1A, JAK1B, Janus kinase 1, JTK3, Tyrosine protein kinase JAK1
Observed MW:	130 kDa
Uniprot ID :	P23458

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:6000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab09903(JAK1 antibody) at dilution of 1:200



Various lysates were subjected to SDS PAGE followed by western blot with FNab09903(JAK1 antibody) at dilution of 1:2000