

anti- CDR2 antibody

Product Information

Catalog No.:	FNab09905
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	cerebellar degeneration-related protein 2
Synonyms:	CDR2, CDR62, PCD17, Yo
Observed MW:	None
Uniprot ID :	Q01850

Application

Reactivity:	Human, rat
Tested Application:	IHC
Recommended dilution:IHC:	1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab09905(CDR2 antibody) at dilution of 1:100